

CAMERA

D-CAM6E



1. Camera

- 4K Resolution
- Auto Focus
- OIS(Optical Image Stabilizer)
- Deeper focus function(25mm)
- View angle : 20° (D-CAM6E), 35° (D-CAM6)
- 2X Digital Zoom
- 0.39" View Finder(FHD)
- Weight : 51g

2. Cable

- From Camera to Transmitter
- Length : 1.5m
- Weight : 53g

3. Transmitter

- Rechargeable Battery 10,000mAh
- Operation time : > 6hrs
- Wireless transfer image to Monitor (HDMI)
- Weight : 157g

4. Foot switch

- Recording / Still shot / 2X Digital zoom

D-CAM4SV



- Wireless HDMI
- 2K / FHD
- Manual Focus
- EIS
- View angle : 30°
- 1.5" View Finder
- Remote Control

D-CAM4S



- Wireless HDMI
- 2K / FHD
- Manual Focus
- EIS
- View angle : 30°
- Remote Control

PRODUCT LINE-UP

PRODUCT	MODEL	SPECIFICATION	PICTURE	PRICE (End User, USD)	REMARKS
HEAD LAMP	DKH-60/60S + DKT-3A or DKT-4A or DKT-6A	- DKT-3A : 2.5X - DKT-4A : 4X - DKT-6A : 6X			Including 2 Batteries 1 Charger set 3 Pairs filter : Transparent, Yellow Brown (for Resin)
	DKH-60/60S + DKT-2 or DKL-6	- DKT-2 : 2X - DKL-6 : 1.5X			
HEAD LAMP	DKH-50/50S	- Illuminance: ~70,000lux - DKH-50S: Motion sensor			
	DKH-40	- Illuminance: 40,000lux			
	DKH-30	- Illuminance: 25,000lux			
LOUPE	DKL-5	- Presbyopia correction loupe - ~1.5X			
	DKT-3	- Galilean Telescope - 2.0X - PD & WD adjustable			
	DKT-4	- Keplerian Telescope - 4.0X - PD & WD adjustable			
	DKT-6	- Keplerian Telescope - 6.0X - PD & WD adjustable			
CAMERA	D-CAM6 D-CAM6E	- 4K Resolution - Auto focus, OIS - View angle: 20 or 35° - View finder - Foot Pedal			
	D-CAM4SV	- 2K / FHD - Manual focus, EIS - View angle: 30° - View finder - Remote control			
	D-CAM4S	- 2K / FHD - Manual focus, EIS - View angle: 30° - Remote control			



WORLD BEST

www.dr-kim.net www.dr-kim.co.kr

ERGONOMIC LOUPE WITH HEADLAMP

• Headlamp(DKH-60/60S)

- Weight : 155g(with a battery)
- Illuminance : ~70,000lux
- Motion sensor : DKH-60S

• Ergonomic Loupes

- DKT-3A: 2.5X, 28g
- DKT-4A: 4X, 32g
- DKT-6A: 6X, 33g

• Direct Loupes(option)

- DKT-2: 2X, 15g
- DKL-6: 1.5X, 8g



2.5X through DKT-3A
4X through DKT-4A
6X through DKT-6A



2X through DKT-2
1.5X through DKL-6

SHADOWLESS HEAD LAMP

DKH-50/50S



- Illuminance: ~70,000lux @ 30cm
- Operation continuity: 2.3~12hours
- Weight: 154g(with a battery)
- Motion sensor: DKH-50S

DKH-40



- Illuminance: 40,000lux
- Operation continuity : 2.5hours
- Weight : 126g (with a battery)

DKH-30



- Illuminance: 25,000lux
- Operation continuity: 5hours
- Weight: 122g(with a battery)

PRESBYOPIA CORRECTION LOUPE

DKL-5



- Magnification: ~1.5X
- For ENT, dentistry or plastic surgery
- Aspherical lenses
- Weight : 16g



ATTACHMENT
Insert the protrusion of the DKL-5 into lower hole of Dr-Kim head lamp

HIGH POWER LOUPE

DKT-3



- Magnification : 2.0X
- Working distance : Adjustable(25~50cm)
- PD : Adjustable(56~72mm)
- Weight : 25g

DKT-6



- Magnification : 6.0X
- Working distance : 40cm, Adjustable
- PD : Adjustable(56~72mm)
- Weight : 37g

DKT-4



- Magnification : 4.0X
- Working distance : 50cm, Adjustable
- PD : Adjustable(56~72mm)
- Weight : 25g



ATTACHMENT
Insert the protrusion of the DKT-3, 4 into upper hole of Dr-Kim head lamp