

CERTIFICATE OF ACCEPTANCE


Serial #	25001529, 25001978 to 25001986, 25002228 to 25002285
Part Number	PC12(39,10,1.25)
Material	Poly(ε-caprolactone)
Date of Expiry	2028-01

Attribute	Acceptance Criteria	Analysis Result
Appearance (pattern design)	Adheres to pattern as described in technical specifications.	Pass
Appearance (Colour)	White to light yellow	Pass
Dimensions	± 10% or 2mm, whichever is smaller	Pass
Sterility Dose for SAL 10 ⁻⁶	16.2-19.8 kGy	Pass

STORAGE INSTRUCTIONS

1. The material remains stable according to its expiry date if kept under the recommended storage conditions.
2. Store in cool dry place out of direct sunlight.
3. Prior to use inspect package for signs of tampering, damage or water contamination.
4. Do not leave the device in the car under direct sunlight even for short periods to prevent accidental melting of the device.

Released by:


Senior Manager
Quality Assurance



Date:

17 MAR 2025

Signature by QA signifies that the batch record has been reviewed and that the above values accurately reflect the results obtained.

Warranty: Disclaimer: We warrant that our products meet all applicable specifications as shown. WE MAKE NO OTHER WARRANTIES, EXPRESS OR IMPLIED, AND SEPCIFICALLY DISCLAIM ALL OTHER WARRANTIES INCLUDING QUALITY, MERCHANTABILITY, FITNESS FOR PURPOSE, OR NON-INFRINGEMENT.