

**LIVE A
VIVACIOUS
LIFE**



**VIVALIV®
VIVARO**

Medical Aesthetic
Medical Device
Customized Cannula



Live a vivacious life

Dongwoo VIVARO


Dongwoo Medical Techno LTD is a specialized medical device manufacturer dedicated to the development and production of high-quality disposable needle products for use in the medical and aesthetic industries.

Since our founding in 2012, we have focused on a wide range of specialized needle - with particular expertise in aesthetic-purpose needles. Given the technical specificity of our products, we continuously invest in research and development to create innovative solutions and respond quickly to the evolving needs of the market. Our flexible organizational structure allows us to work closely and responsively with our clients to deliver tailored solutions.

VIVARO Cannula

Vivaro Cannula is a high-precision, single-use blunt cannula designed for aesthetic procedures involving dermal and subdermal injections. It enables smooth and controlled delivery of high-viscosity substances such as hyaluronic acid fillers while minimizing patient discomfort, bruising, and tissue trauma. Manufactured using DONGWOO M TECHNO's advanced automated machining technology, Vivaro Cannula ensures consistent outer diameter quality, smooth insertion, and exceptional handling performance, allowing practitioners greater precision and control during treatment.

Ideal for procedures including hyaluronic acid filler injection, fat grafting, and mesotherapy, Vivaro Cannula enhances patient satisfaction by reducing pain, swelling, and recovery time, especially in sensitive areas such as under the eyes, lips, and nasolabial folds.



Gauge	Length(mm)	Gauge	Length(mm)	Gauge	Length(mm)
16G	100 / 150	20G	40 / 50 / 60 / 70 / 90 / 100 / 120 / 150	24G	40 / 50 / 70
17G	50 / 70 / 90 / 100 / 150	21G	40 / 50 / 60 / 70 / 90 / 100	25G	30 / 40 / 50 / 70
18G	50 / 70 / 90 / 100 / 150	22G	40 / 50 / 70	27G	30 / 40 / 50 / 70
19G	40 / 50 / 60 / 70 / 90 / 100	23G	40 / 50 / 70 / 80	30G	25 / 30 / 38



VIVAR Cannula


Model No.	SIZE
DUB01	16G/100mm
DUB06	17G/100mm
DUB12	18G/100mm
DUB37	21G/50mm
DUB39	21G/70mm
DUB43	22G/50mm
DUB44	22G/70mm
DUB56	23G/50mm
DUB57	23G/70mm
DUB71	25G/40mm
DUB72	25G/50mm
DUB73	25G/70mm
DUB75	27G/40mm
DUB76	27G/50mm
DUB85	30G/25mm

Customized cannula

At DONGWOO M TECHNO, we understand that each medical application demands a unique solution. That's why we specialize in the custom manufacturing of precision needles, meticulously designed to meet the specific requirements of our clients.

From gauge size and length to bevel angles and hub compatibility, every element is engineered with accuracy and flexibility. Whether you're developing a new medical device or require optimization for an existing one, our advanced production capabilities and quality-driven process ensure that your vision becomes a reality with unmatched precision.

12 hole Irrigation cannula



Customized cannula types

D-Type



L-Type




W-Type



WA-Type



Raw Material Info



Posco
Stainless Steel 304



Dupont
tyvek 1073B CR27 / 2FS



Pet
Medical Blister Package

Production Manufacturing Information



Clean room
ISO 14644-1
Class 8 or 100,000 : 0.5µm/ft³



Sterile Pouch

Adress

208, 209, Techno zone,
117 Hwanggeum-ro, Yangchon-Eup,
Gimpo-si, Gyeonggi-do, Korea

Contact Information

TEL +82 31 981 3448
FAX +82 31 981 3449

Customer service

jihoonlee@dmcare.co.kr
+82 10 8227 8408

www.dongwoomtechno.com

